

Thermo Scientific DreamTaq DNA Polymerase

Catalog Number EP0701, EP0702, EP0703, EP0704, EP0705

Pub. No. MAN0012036 Rev. C.00



WARNING! Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Sheets (SDSs) are available from thermofisher.com/support.

Contents and storage

Cat. No.	Contents	Amount	Storage
EP0701	DreamTaq DNA Polymerase, 5 U/μL	200 U	-25 °C to -15 °C
	10X DreamTaq Buffer*	1.25 mL	
EP0702	DreamTaq DNA Polymerase, 5 U/μL	500 U	
	10X DreamTaq Buffer*	2 x 1.25 mL	
EP0703	DreamTaq DNA Polymerase, 5 U/μL	5 x 500 U	
	10X DreamTaq Buffer*	10 x 1.25 mL	
EP0704	DreamTaq DNA Polymerase, 5 U/μL	20 x 500 U	
	10X DreamTaq Buffer*	40 x 1.25 mL	
EP0705	DreamTaq DNA Polymerase, 5 U/μL	10 x 500 U	
	10X DreamTaq Buffer*	20 x 1.25 mL	

* includes 20 mM MgCl₂

Description

Thermo Scientific™ DreamTaq™ DNA Polymerase is an enhanced *Taq* DNA polymerase optimized for all standard PCR applications. It ensures higher sensitivity, longer PCR products and higher yields compared to conventional *Taq* DNA polymerase.

DreamTaq DNA Polymerase uses the same reaction set-up and cycling conditions as conventional *Taq* DNA polymerase. Extensive optimization of reaction conditions is not required. The enzyme is supplied with optimized DreamTaq buffer, which includes 20 mM MgCl₂.

DreamTaq DNA Polymerase generates PCR products with 3'-dA overhangs. The enzyme is inhibited by dUTP but can incorporate modified nucleotides.

Features

- Robust amplification with minimal optimization.
- High yields of PCR products.
- Higher sensitivity compared to conventional *Taq* DNA polymerase.
- Amplification of long targets up to 6 kb from genomic DNA and up to 20 kb from viral DNA.
- Generates 3'-dA overhangs.
- Incorporates modified nucleotides.

Applications

- Routine PCR amplification of DNA fragments up to 6 kb from genomic DNA and up to 20 kb from viral DNA.
- RT-PCR.
- Generation of PCR products for TA cloning.

Definition of Activity Unit

One unit of the enzyme catalyzes the incorporation of 10 nmol of deoxyribonucleotides into a polynucleotide fraction in 30 min at 74 °C.

10X DreamTaq Buffer

DreamTaq Buffer is a proprietary formulation which contains KCl and (NH₄)₂SO₄ at a ratio optimized for robust performance of DreamTaq DNA Polymerase in PCR applications. DreamTaq Buffer also includes MgCl₂ at a concentration of 20 mM.

For Research Use Only. Not for use in diagnostic procedures.

ThermoFisher
SCIENTIFIC

Inhibition and Inactivation

- Inhibitors: ionic detergents (deoxycholate, sarkosyl and SDS) at concentrations higher than 0.06, 0.02 and 0.01 %, respectively.
- Inactivated by phenol/chloroform extraction.

Protocol

To set up parallel reactions and to minimize the possibility of pipetting errors, prepare a PCR master mix by mixing water, buffer, dNTPs, primers and DreamTaq DNA Polymerase. Prepare enough master mix for the number of reactions plus one extra. Aliquot the master mix into individual PCR tubes and then add template DNA.

1. Gently vortex and briefly centrifuge all solutions after thawing.
2. Place a thin-walled PCR tube on ice and add the following components for each 50 µL reaction:

Contents	Volume
10X DreamTaq Buffer*	5 µL
dNTP Mix, 2 mM each (#R0241)	5 µL (0.2 mM of each)
Forward primer	0.1-1.0 µM
Reverse primer	0.1-1.0 µM
Template DNA	10 pg - 1 µg
DreamTaq DNA Polymerase	1.25 U
Water, nuclease-free (#R0581)	to 50 µL
Total volume	50 µL

*10X DreamTaq Buffer contains 20 mM MgCl₂, which is optimal for most applications. If additional optimization is required, 25 mM MgCl₂ (#R0971) can be added to the master mix. The volume of water should be reduced accordingly.

Volumes of 25 mM MgCl₂, required for specific final MgCl₂ concentration:

Final concentration of MgCl ₂	2 mM	2.5 mM	3 mM	4 mM
Volume of 25 mM MgCl ₂ to be added for 50 µL reaction	0 µL	1 µL	2 µL	4 µL

3. Gently vortex the samples and spin down.
4. When using a thermal cycler that does not contain a heated lid, overlay the reaction mixture with 25 µL of mineral oil.
5. Place the reactions in a thermal cycler. Perform PCR using the recommended thermal cycling conditions outlined below:

Step	Temperature, °C	Time	Number of cycles
Initial denaturation	95	1-3 min	1
Denaturation	95	30 s	25-40
Annealing	T _m -5	30 s	
Extension*	72	1 min	
Final Extension	72	5-15 min	1

* The recommended extension step is 1 min for PCR products up to 2 kb. For longer products, the extension time should be prolonged by 1 min/kb.

Guidelines for Preventing Contamination of PCR Reaction

During PCR more than 10 million copies of template DNA are generated. Therefore, care must be taken to avoid contamination with other templates and amplicons that may be present in the laboratory environment. General recommendations to lower the risk of contamination are as follows:

- Prepare your DNA sample, set up the PCR mixture, perform thermal cycling and analyze PCR products in separate areas.
- Set up PCR mixtures in a laminar flow cabinet equipped with an UV lamp.
- Wear fresh gloves for DNA purification and reaction set up.
- Use reagent containers dedicated for PCR. Use positive displacement pipettes or use pipette tips with aerosol filters to prepare DNA samples and perform PCR set up.
- Use PCR-certified reagents, including high quality water (e.g., Water, nuclease-free, #R0581).
- Always perform "no template control" (NTC) reactions to check for contamination.

DreamTaq DNA Polymerase does not incorporate dUTP, therefore it is not possible to perform carryover contamination prevention with UDG. For this application we recommend using Taq DNA Polymerase (#EP0401) or Thermo Scientific™ DreamTaq™ Hot Start DNA Polymerase (#EP1701).

Guidelines for Primer Design

Use primer design software or follow the general recommendations for PCR primer design as outlined below:

- PCR primers are generally 15-30 nucleotides long.
- Optimal GC content of the primer is 40-60 %. Ideally, C and G nucleotides should be distributed uniformly along the primer.
- Avoid placing more than three G or C nucleotides at the 3'-end to lower the risk of non-specific priming.
- If possible, the primer should terminate with a G or C at the 3'-end.
- Avoid self-complementary primer regions, complementarities between the primers and direct primer repeats to prevent hairpin formation and primer dimerization.
- Check for possible sites of undesired complementarity between primers and template DNA.
- When designing degenerate primers, place at least 3 conserved nucleotides at the 3'-end.
- When introducing restriction enzyme sites into primers, refer to the table "Reaction conditions for FastDigest enzymes" located on www.thermofisher.com/fastdigest to determine the number of extra bases required for efficient cleavage.
- Differences in melting temperatures (T_m) between the two primers should not exceed 5 °C.

Estimation of Primer Melting Temperature

For primers containing less than 25 nucleotides, the approx. melting temperature (T_m) can be calculated using the following equation:

$$T_m = 4 (G + C) + 2 (A + T),$$

where G, C, A, T represent the number of respective nucleotides in the primer.

If the primer contains more than 25 nucleotides specialized computer programs are recommended to account for interactions of adjacent bases, effect of salt concentration, etc.

Components of the Reaction Mixture

Template DNA

Optimal amounts of template DNA for a 50 µL reaction volume are 0.01-1 ng for both plasmid and phage DNA, and 0.1-1 µg for genomic DNA. Higher amounts of template increase the risk of generation of non-specific PCR products. Lower amounts of template reduce the accuracy of the amplification.

Trace amounts of certain agents used for DNA purification, such as phenol, EDTA and proteinase K, can inhibit DNA polymerases. Ethanol precipitation and repeated washes of the DNA pellet with 70% ethanol normally removes trace contaminants from DNA samples.

MgCl₂ Concentration

DreamTaq DNA Polymerase is provided with an optimized 10X DreamTaq Buffer which includes MgCl₂ at a concentration of 20 mM. A final MgCl₂ concentration of 2 mM is generally ideal for PCR. The MgCl₂ concentration can be further increased up to 4 mM by the addition of 25 mM MgCl₂ (#R0971).

If the DNA samples contain EDTA or other metal chelators, the Mg²⁺ ion concentration in the PCR mixture should be increased accordingly (1 molecule of EDTA binds 1 Mg²⁺).

dNTPs

The recommended final concentration of each dNTP is 0.2 mM. In certain PCR applications, higher dNTP concentrations may be necessary. It is essential to have equal concentrations of all four nucleotides (dATP, dCTP, dGTP and dTTP) present in the reaction mixture.

To obtain a 0.2 mM concentration of each dNTP in the PCR mixture, please refer to the table below:

Volume of PCR mixture	dNTP Mix, 2 mM each (#R0241)	dNTP Mix, 10 mM each (#R0191)	dNTP Mix, 25 mM each (#R1121)
50 µL	5 µL	1 µL	0.4 µL
25 µL	2.5 µL	0.5 µL	0.2 µL
20 µL	2 µL	0.4 µL	0.16 µL

Primers

The recommended concentration range of the PCR primers is 0.1-1 µM. Excessive primer concentrations increase the probability of mispriming and generation of non-specific PCR products.

For degenerate primers and primers used for long PCR, we recommend higher primer concentrations in the range of 0.3-1 µM.

Cycling Parameters

Initial DNA Denaturation

It is essential to completely denature the template DNA at the beginning of the PCR run to ensure efficient utilization of the template during the first amplification cycle. If the GC content of the template is 50 % or less, an initial 1-3 min denaturation at 95 °C is sufficient. For GC-rich templates this step should be prolonged up to 10 min. If a longer initial denaturation step is required, or if the DNA is denatured at a higher temperature, DreamTaq DNA Polymerase should be added after the initial denaturation step to avoid a decrease in its activity.

Denaturation

A DNA denaturation time of 30 seconds per cycle at 95 °C is normally sufficient. For GC-rich DNA templates, this step can be prolonged to 3-4 min. DNA denaturation can also be enhanced by the addition of either 10-15 % glycerol, 10 % DMSO, 5 % formamide or 1-1.5 M betaine. The melting temperature of the primer-template complex decreases significantly in the presence of these reagents. Therefore, the annealing temperature has to be adjusted accordingly.

In addition, 10 % DMSO and 5 % formamide inhibit DNA polymerases by 50 %. Thus, the amount of the enzyme in the reaction should be increased if these additives are used.

Primer Annealing

The annealing temperature should be 5 °C lower than the melting temperature (T_m) of the primers. Annealing for 30 seconds is normally sufficient. If non-specific PCR products appear, the annealing temperature should be optimized stepwise in 1-2 °C increments. When additives which change the melting temperature of the primer-template complex are used (glycerol, DMSO, formamide and betaine), the annealing temperature must also be adjusted.

Extension

The optimal extension temperature for DreamTaq DNA Polymerase is 70-75 °C. The recommended extension step is 1 min at 72 °C for PCR products up to 2 kb. For longer products, the extension time should be prolonged by 1 min/kb. For amplification of templates >6 kb a reduction of the extension temperature to 68 °C is recommended to avoid enzyme inactivation during prolonged extension times.

Number of Cycles

The number of cycles may vary depending on the amount of template DNA in the PCR mixture and the expected PCR product yield. If less than 10 copies of the template are present in the reaction, about 40 cycles are required. For higher template amounts, 25-35 cycles are sufficient.

Final Extension

After the last cycle, it is recommended to incubate the PCR mixture at 72 °C for additional 5-15 min to fill-in any possible incomplete reaction products. If the PCR product will be cloned into TA vectors, the final extension step may be prolonged to 30 min to ensure the complete 3'-dA tailing of the PCR product. If the PCR product will be used for cloning using Thermo Scientific™ CloneJET™ PCR Cloning Kit (#K1231), the final extension step can be omitted.

Troubleshooting

For troubleshooting, please visit www.thermofisher.com

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale at www.thermofisher.com/us/en/home/global/terms-and-conditions.html. If you have any questions, please contact Life Technologies at www.thermofisher.com/support.



Thermo Fisher Scientific Baltics UAB | V.A. Graiciuno 8, LT-02241 Vilnius, Lithuania

For descriptions of symbols on product labels or product documents, go to thermofisher.com/symbols-definition.

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2022 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

thermofisher.com/support | thermofisher.com/askaquestion

thermofisher.com

ThermoFisher
S C I E N T I F I C

28 April 2022

